




™ Simple
Connected
Powerful



MTS101 Monitor with Touch Screen

The MTS101 is a medical 10.1" USB powered Multi-Touch Monitor that can be used in the OR for controlling all MediCapture MVR imaging systems.

Designed for modern operating rooms and surgery centers, the MTS101 enhances the user experience by providing perfect monitor positioning and a larger user interface for entering patient data, editing videos and developing reports.



Designed for the OR

The MTS101 eliminates operational boundaries. You can now control the MVR system in any optimal location within the OR. You will also benefit from a larger user interface to capture, edit, report and share video footage instantaneously.


Easy fit

With its sizable 10.1" (26cm) touch screen and easy VESA Mount capability, you can position the MTS101 in any convenient location throughout the OR.


Simple to Use

With MTS101 and its 5-points Multi Touch, you can handle your workflow with a swipe of finger -- just like with your smartphone.

MTS101 Features



Backside with VESA 75



Side connections

Medical monitor with touch screen and fully sealed metal housing.

Slim design, solid screen protection.

Plug and play setup. Connect the monitor to the USB and HDMI, and it will operate when your host device is turned on.

Multiple ergonomic mounting options with VESA 75.

MTS101 Specs¹

Device Type	Medical Monitor with Touch Screen Aluminium Alloy, 26cm / 10.1" diagonal, 260 cd/m ² , Aspect Ratio: 16:10, Multi-Touch 5 points
Monitor Resolution	1920 x 1200
Video Inputs	MINI HDMI (type C) Supports HDMI (DVI via adapter)
Input Resolution	1920x1200p60
USB Support	USB type C 2.0
Mechanical	Size: 244,1 x 167,1 x 20,2 mm / 9,6 x 6,6 x 0,8" Weight: 0,63 kg (1,4lbs)
Environment	Storage: -40° to +85° C Operating: -20°C to +40°C
Power	From USB-C, DC 5V / 0.8A
Certifications	EN 60601-1:2006 + A1:2013 + A12:2014 - Medical Equipment Safety, UL, MDR 2017/745 EU, Class 1, EN 60601-1-2:2015 - Electromagnetic Compatibility EMC, FCC Part 15B - Class B, RoHS, WEEE, UKCA
Order Number	MTS101 : MTS monitor, HDMI to MINI HDMI cable, USB A to USB C cable, Information for Use (English) ² ST004MTS : Stand, Acryl, for MTS101



¹ Subject to change without notice

² For more languages of the User Guide visit www.medicapture.com/support

MediCapture, Inc. - North America

2250 Hickory Road, Suite 200
Plymouth Meeting, PA 19462, USA

+1 888 922 7887

www.medicapture.com
sales@medicapture.com

MediCapture Europe B.V.

Kraanpoort 10
6041 EG Roermond
The Netherlands

+31 475 214 305

www.medicapture.com
sales.europe@medicapture.com

MediCapture APAC Inc.

3FL-6, No. 432, Section 1, Keelung Road,
Xinyi District 110, Taipei City,
Taiwan, Republic Of China

+886 2 2345 3668

www.medicapture.com
sales.asiapacific@medicapture.com