



™ Simple
Connected
Powerful



MVR Pro

Powered by the exclusive iMave Pro platform, the MVR Pro is our most advanced, premium-quality medical video recorder. With full-patient workflow capability and an easy built-in touchscreen, this advanced recorder is designed to benefit today's modern operating rooms and surgery centers.



The Next-Generation Video Recorder

Multi-featured, the MVR Pro benefits from MediCapture's unique Smart Workflow that enables surgical teams to control the recorder according to preference including touchscreen or tablets⁴. It also enables the user to jump from recording to managing data with the swipe of a finger.

Full Patient Workflow

The MVR Pro comes with full patient workflow, including a DICOM Upgrade Package¹ option, for easy and fast retrieval of patient information and transfer of studies to PACS.

Simple to Use

The unique built-in touchscreen includes intuitive icons making the recorder very easy to use. No formal training is required.

Unique Features



MVR Manager

Browse and review all patient files. Zoom, take snapshots, cut video. Add notes, report in PDF, print and archive as you wish.



MVR Secure

Protect your files from hacking or power outages. Includes user authentication and historical logs for data protection.

Upgrade Package

Install MVR Upgrade Packages at any time by simply entering the activation codes provided by MediCapture.



DICOM¹

Easily retrieve patient DICOM worklist and store to your PACS.



4K Recording^{1,5}

Records delicate surgeries in stunning 4K (2160 p30/25).

Also available : MVR

MVR offers same features set as MVR Pro without the integrated screen. It can be operated from external USB touchscreen, Android tablet⁴ or keypad/mouse.



MVR Remote App⁴



- ¹ Upgrade package option
- ² Subject to changes without notice
- ³ Please ask for a quotation for higher memory
- ⁴ Install the MVR Remote App on your Android smartphone or tablet
- ⁵ Check camera compatibility with MediCapture

MVR Pro / MVR Specs²

Device Type	HD Medical Video Recorder Records High-Definition video and still images from medical video sources to: USB, internal storage, network storage. Exports video and images to PACS in DICOM format ¹
Case/Buttons	Stainless steel, fluid resistant
User Interface	Overlay on internal touchscreen monitor, or external monitor
Video Inputs	Multi-Input: 2x HDMI, DVI (with adapter), 4K-SDI with MVC Pro SDI to HDMI (optional) ⁵
Video Outputs	1x HDMI, DVI (with adapter)
Video Input Resolutions	4K on Input 1 only ⁵ : 3840x2160p30 Both inputs (1+2): 1920x1200p60, 1080p60/50/30/25, 1080i60/50, 1280x1024p60/50, 720p 60/50
Video Recording Resolutions	3840x2160 ¹ , 1920x1200, 1920x1080, 1280x1024, 1280x720. For 3840x2160 recording, "4K recording" needs to be activated ^{1,5}
Video Recording Bitrates	MP4 / coding H.264 up to 32 Mbit/s (1920x1080 p60), lower resolutions have lower bitrates MP4 / coding HEVC (H.265) reduces bitrate by 50%.
Image Format	JPG, PNG
Sound Recording	Bluetooth Mic, Mic by Stereo Headset (4 poles), USB Mic / Audio embedded in video
Internal Storage	128GB. Extension on request ³
External Storage	USB flash drive, USB hard drive (FAT32, NTFS, exFAT)
USB Support	USB3.0, SuperSpeed, 1 front connector, 2 rear connectors
Network Storage	External Server, 1000/100/10Mbit/s
Footswitch/Trigger	3,5mm mini jack. Compatible with standard footswitches (not included) and closed-contact camera triggers, Bluetooth footswitches
Remote Control	IR remote (included), RS232 (via USB adapter), JSON API over http and over RS232, MVR Remote App ⁴
Mechanical (touchscreen closed)	Size: 173 x 153 x 50mm / 6,7 x 6,0 x 2" Weight: 1,20 kg / 2lb 13oz
Environmental	Storage: -40°C to +85°C Operating: -20°C to +40°C
Power	AC 100-240V, 50-60Hz, 0.4-0.2A, 3 prong AC jack
Certifications	MDR 2017/745, CE 93/42/EEC, EN 60601-1:2006 + A1:2013 + A12:2014 (Ed 3.1), EN 60601-1-2:2015 (Edition 4), EN 60601-1-6:2010+A1:2015, EN 50581:2012, RED 2014/53/EU, FCC Part 15, Subpart B, Class B Class 1 Equipment according to US-FDA and MDR 2017/745, RoHS3 / WEEE / ISO 13485: 2016
	For more languages of the User Guide visit www.medicapture.com/support
Order Number	MVR409 / MVR420 MVR recorder, USB flash drive, HDMI cable, Wi-Fi / BT antenna, HDMI to DVI adapter, Split Remote Cable Adapter, Power Cable, IR Remote Control, Rubber Feet (4), User Guide (English)
Order information for optional Upgrade Packages:	AKDICOM : DICOM Upgrade Package AK4K : 4K Recording Upgrade Package

MediCapture, Inc. – North America

2250 Hickory Road, Suite 200
Plymouth Meeting, PA 19462, USA

+1 888 922 7887

www.medicapture.com
sales@medicapture.com

MediCapture Europe B.V.

Stationsplein 5
6131 AT Sittard
The Netherlands

+31 46 208 1120

www.medicapture.com
sales.europe@medicapture.com

MediCapture APAC Inc.

3FL-6, No. 432, Section 1, Keelung Road,
Xinyi District 110, Taipei City,
Taiwan, Republic Of China

+886 2 2345 3668

www.medicapture.com
sales.asiapacific@medicapture.com