

4K p60



™ Simple
Connected
Powerful




MVR460

Medical Video Recorder
4K


USB (PTZ) Camera Interface¹



USB Endoscope Interface¹



Universal Printer Compatibility²




Archive & Send to PACS/EHR



Access to Patient Worklist




Streaming & Sharing



- DICOM
- 12G SDI
- DUAL INPUT
- SMART WORKFLOW
- 2TB
- MULTI STORAGE
- BUILT-IN SPEAKERS

Medical Imaging Made Easy™

Use Cases



Medical Imaging Environment (OR)

Makes it easy to record and archive medical imaging and manage the Patient Worklist.



Robot-Assisted Surgical Units

Ensures maximum efficiency in robot-assisted surgeries with an easy-mount design for quick installation.

MVR460 Specifications⁴

| | |
|---|---|
| Device Type | Medical Video Recorder |
| Case & Buttons | Aluminium Alloy Case, VESA100 on bottom, Stainless Steel Buttons (Fluid Resistant) |
| Video Connections (Accepted Input signals) | 2 selections of: HDMI 4K, DisplayPort 4K, USB camera 4K, SDI up to 12G-SDI |
| Recording Resolution | 3840x2160, 1920x1080, 1280x1024, 1280x720, 1024x768 |
| Video Recording Format | MP4: H.264 / HEVC H.265 High Efficiency Video Coding |
| Image format | JPG / PNG |
| Audio Input (Sound recording) | MIC-in (low level), 3.5mm Mono TS (2-pin), USB Microphone. Line-in 1Vpp, 3.5mm Stereo TRS (3-pin), HDMI |
| Audio Output | Integrated Speaker. LINE-out 1Vpp, 3.5mm Stereo TRS (3-pin). Headset Output: 3.5mm Stereo TRS (3-pin) |
| Internal Storage | HDD 2TB (Larger internal storage by order) |
| External Storage | USB drives (flash drives or hard drives), FAT32 / NTFS / exFAT |
| USB Support | USB3.0, Super speed, 1 front connector, 2 rear connectors |
| Remote | 3.5mm mini jack. Compatible with standard footswitches (not included) and closed-contact camera triggers |
| LAN | RJ-45, Ethernet, max 2.5 GB/s |
| Mechanical (w x d x h) | Dimensions w x h x d (monitor closed): 195 x 71 x 220.3 mm (7,7 x 2,8 x 8,7") |
| Monitor dimensions | 7" diagonal |
| Weight | MVR460: 2,1 kg (4.63 lbs) |
| Ingress protection | IP40 |
| Power | AC 100-240V, 0.6-0.3A, 50/60Hz |
| Box Contents | MVR460 recorder, USB flash drive, HDMI cable, DisplayPort cable, 12G-SDI cable (BNC), Split Remote Cable Adapter, Audio cable M-M, Audio cable M-F, Network Patch Cable RJ45, Power cable, Rubber Feet (4x1) (2), Information for Use (English ⁶) |
| Compliance and Certification | EN 60601-1 - Medical Equipment Safety EN 60601-1-2 - EMC, EN IEC 63000:2018 - Restriction of hazardous substances CE (MDR 2017/745/EU), UKCA, WEEE, RoHS, FCC, GDPR and HIPAA compliant. |
| Classification | Class I Equipment / Mode of operation: continuous No AP/APG. No Applied Part |
| Order Information | MVR460 Medical Video Recorder AKDICOM - DICOM Software activation |

MediCapture MVR REMOTE APP



Optional add-ons Pay for what you really need



¹ Please check with MediCapture for known compatible cameras and devices.
² This device can connect with most network-based printers.
³ Upgrade package option.
⁴ These specifications are subject to change without notice.
⁵ Will be available in the future.
⁶ For more languages visit www.medicapture.com/support

MediCapture, Inc. - North America

2250 Hickory Road, Suite 200
Plymouth Meeting, PA 19462, USA

+1 888 922 7887
www.medicapture.com
sales@medicapture.com

MediCapture Europe B.V.

Stationsplein 5 C2
6131 AT Sittard, The Netherlands

+31 46 208 1120
www.medicapture.com
sales.europe@medicapture.com

MediCapture APAC Inc.

3FL-6, No. 432, Section 1, Keelung Road,
Xinyi District 110, Taipei City, Taiwan, Republic Of China

+886 2 2345 3668
www.medicapture.com
sales.asiapacific@medicapture.com